


# Next Generation Sequencing (NGS)

HAMIDAH SUUBI NAMAGEMBE

26<sup>TH</sup> JUNE 2023


# Introduction

- Next-Generation Sequencing (NGS) is a high-throughput DNA sequencing method that has transformed genomic research
- NGS enables faster, more cost-effective, and more comprehensive analysis of genetic material




# NGS Workflow

## STEP 1: Extraction



## STEP 2: Library Prep



## STEP 3: Sequencing



## STEP 4: Analysis



# NGS Technologies

## Short-Read NGS

- Sequencing by synthesis (SBS) –
  - i. Cyclic reversible chain termination (CRT):* Illumina & Qiagen
  - ii. Single nucleotide addition (SNA):* Roche & Ion Torrent
- Sequencing by ligation (SBL)- Thermo Fisher & BGI


## Long-Read NGS

- Realtime sequencing approaches – Pacific Biosciences & Oxford Nanopore Technologies
- Synthetic approaches- Illumina & 10X Genomics

# Short-Read NGS


Sequencing by Synthesis (SBS)- *Cyclic reversible chain termination* (CRT): Illumina

a Illumina




## Nucleotide addition

Fluorophore-labelled, terminally blocked nucleotides hybridize to complementary base. Each cluster on a slide can incorporate a different base.



## Imaging

Slides are imaged with either two or four laser channels. Each cluster emits a colour corresponding to the base incorporated during this cycle.




## Cleavage

Fluorophores are cleaved and washed from flow cells and the 3'-OH group is regenerated. A new cycle begins with the addition of new nucleotides.

# Short-Read NGS


## Sequencing by ligation (SBL)- Thermo Fisher

**a SOLiD**  
(Thermo Fisher)




### Two-base-encoded probes

Probes with two known bases followed by degenerate or universal bases hybridize to a template; ligase immobilizes the complex and the slide is imaged




### Cleavage

The fluorophore is cleaved from the probe along with several bases, revealing a 5' phosphate



### Probe extension

10 rounds of hybridization, ligation, imaging and cleavage identify 2 out of every 5 bases



### Reset


After a round of probe extension, all probes and anchors are removed and the cycle begins again with an offset anchor



# Long-Read NGS

## Single-molecule sequencing approaches-

**ONT** is based on the detection of changes in ionic current as DNA or RNA molecules pass through a nanopore




# Long-Read NGS

**Synthetic approaches-10X  
Genomics -Gel Bead-in-Emulsion  
(GEM) Generation:** DNA  
molecules are encapsulated  
within droplets

## Bb 10X Genomics

### Emulsion PCR

Arbitrarily long DNA  
is mixed with beads  
loaded with  
barcoded primers,  
enzyme and dNTPs




### GEMs

Each micelle  
has 1 barcode  
out of 750,000




### Amplification

Long fragments are  
amplified such that the  
product is a barcoded  
fragment ~350 bp



### Pooling

The emulsion is  
broken and DNA is  
pooled, then it  
undergoes a standard  
library preparation



### Linked reads

- All reads from the same GEM derive from the long fragment, thus they are linked
- Reads are dispersed across the long fragment and no GEM achieves full coverage of a fragment
- Stacking of linked reads from the same loci achieves continuous coverage



# Short-read vs Long-read

**Short-read** NGS involves sequencing DNA fragments of relatively small lengths, typically ranging from 50 to 300 base pairs.


High-throughput: a large volume of sequences in a single run

Cost-effective: generally more affordable per base

High accuracy: Short reads typically have low error rates

- **Long-read NGS**-Long DNA fragments of thousands of base pairs
- Greater sequence continuity: better characterization of complex genomic regions, repetitive sequences, and structural variations
- Direct detection of modifications: such as DNA methylation

# Applications of Short-read and Long-read platforms



- **Long-read**
- De novo genome assembly: for reconstructing complex genomes with high accuracy and resolving structural variations
- Characterizing repetitive regions: help resolve repetitive regions
- Studying RNA isoforms: provides insights into alternative splicing events and full-length isoforms in transcriptomics research

# Data Quality Metrics

- **Depth of coverage:** how many sequence reads are present at a given position
- **Base quality:** have the correct bases been called in sequence reads)
- **Mapping quality:** have the reads been mapped to the correct position in the genome)

Table 1: Quality Scores and Base Calling Accuracy

Phred Quality Score	Probability of Incorrect Base Call	Base Call Accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1,000	99.9%
40	1 in 10,000	99.99%
50	1 in 100,000	99.999%

# Challenges and limitations



**Data Storage:** NGS generates large volumes of data that require substantial storage capacity.



**Computational Requirements:** Complex data analysis necessitates powerful computing



**Cost:** Although prices have decreased, NGS can still be relatively expensive.



**Sample Quality:** Poor sample quality may impact sequencing results and data interpretation.

# References

- doi: [10.1177/1099800417750746](https://doi.org/10.1177/1099800417750746)
- <https://irepertoire.com/ngs-overview-from-sample-to-sequencer-to-results/>
- Goodwin, S., McPherson, J. D., & McCombie, W. R. (2016). Coming of age: ten years of next-generation sequencing technologies. *Nature Reviews Genetics*, 17(6), 333–351. doi:10.1038/nrg.2016.49





THANK YOU

# NGS Technologies

---

## **Short- Read NGS**

Sequencing by synthesis (SBS) –*Cyclic reversible chain termination* (CRT): Illumina & Qiagen, and *Single nucleotide addition* (SNA): Roche & Ion Torrent

---

Sequencing by ligation (SBL)- Thermo Fisher, BGI

---

## **Long- Read NGS**

Realtime sequencing approaches – Pacific Biosciences, Oxford Nanopore Technologies


---

Synthetic approaches- Illumina, 10X Genomics


---

# Short-Read NGS

## Sequencing by Synthesis (SBS)




## Sequencing by Ligation (SBL)



# Long-Read NGS

## Single-molecule sequencing approaches

### Ab Oxford Nanopore Technologies




## Synthetic approaches

### Bb 10X Genomics

#### Emulsion PCR

Arbitrarily long DNA is mixed with beads loaded with barcoded primers, enzyme and dNTPs



#### GEMs

Each micelle has 1 barcode out of 750,000

#### Amplification

Long fragments are amplified such that the product is a barcoded fragment ~350 bp

#### Pooling

The emulsion is broken and DNA is pooled, then it undergoes a standard library preparation



#### Linked reads

- All reads from the same GEM derive from the long fragment, thus they are linked
- Reads are dispersed across the long fragment and no GEM achieves full coverage of a fragment
- Stacking of linked reads from the same loci achieves continuous coverage

# Short-read vs Long-read



**Short-read** NGS involves sequencing DNA fragments of relatively small lengths, typically ranging from 50 to 300 base pairs.



High-throughput: a large volume of sequences in a single run



Cost-effective: generally more affordable per base



High accuracy: Short reads typically have low error rates

- **Long-read NGS**-Long DNA fragments of thousands of base pairs.
- Greater sequence continuity: better characterization of complex genomic regions, repetitive sequences, and structural variations.
- Detection of long-range genomic information: includes intergenic regions, long-range interactions, and haplotype phasing.
- Direct detection of modifications: ONT platforms can detect DNA modifications, such as DNA methylation, directly.